


PROFESSIONAL SOLUTIONS FROM ARDO MEDICAL

SYSTEM FOR URINE COLLECTION AND SPECIMEN SAMPLING

PVC-free Urine bags for the collection of urine and specimen sampling for laboratory testing. Designed with comfort and protection for delicate skin in mind. Our sophisticated urine collection system can be applied to infants (including premature babies) and children of all ages and sexes.

The urine bags can be easily closed with the butterfly shaped adhesive strip, ensuring the urine sample is always ready for the laboratory.



URINE BAGS FOR NEWBORNS/BABIES AND CHILDREN

User friendly

- ✓ Simple to use, time-saving and hygienic collection of a urine sample
- ✓ Secure closure by means of butterfly shaped adhesive area
- ✓ Urine specimen is immediately ready for the laboratory
- ✓ Different shapes and sizes

Safe and hygienic for premature babies, infants and children

- Dual chamber construction with return valve:
 - ✓ Keeps the urine away from the skin
 - ✓ Prevents skin irritation
- High level of skin protection and skin care through:
 - ✓ Perforated and breathable carrier material
 - ✓ Adhesive surface can be easily and safely removed from the skin
 - ✓ Particularly gentle on the skin despite its high adhesive power
 - ✓ Latex-free adhesive surfaces

Cost efficient solution for Hospital use


- The Urine bag allows for a quick collection of a usable urine sample
 - ✓ Unisex design halves the quantity of stock required
 - ✓ Economical use of bags due to good adhesive adhesion
 - ✓ Easy, time-saving and safe application
 - ✓ Low cost per application

Paediatric

Newborn



Preemie



SWISS QUALITY 

Ardo medical Ltd. T +44 (0)1823 336362

www.ardomedical.co.uk

 **ARDO**
Hospital Recommended