

PHILIPS



Oxygen Concentrator – Oxygenate 5

Oxygenate 5 delivers quality care with the patient in mind

Oxygenate 5 is designed to deliver quality patient care and offer an enhanced experience for patients in need of supplemental oxygen.



Delivers a quiet experience with a sound level of <40dB, which is similar to the hum of a refrigerator.^{1,2}



Ships with four additional filters and one additional flat sponge air inlet filter that may help lower consumable purchases and operational cost.



Encourages patients to feel more engaged in their therapy with features like an easy-to-read digital display.



Helps homecare providers and clinicians manage patient care and is supported by customer service that customers can rely on.

¹ Centers for Disease Control and Prevention. (n.d.). Hearing loss and how to prevent it. National Center for Environmental Health. U.S. Department of Health and Human Services. https://www.cdc.gov/nceh/hearing_loss/infographic/

² McKee, G. (2003). Hearing loss in children: Causes and prevention. PubMed Central. <https://pmc.ncbi.nlm.nih.gov/articles/PMC1253729/>



Flowmeter provides a digital reading of oxygen flow that is designed to make reading and setting up parameters easier



User-replaceable filters may help with ease of maintenance



Oxygenate 5 is designed with nine levels of filtration*

Product specifications

Power supply	230V/50Hz
Liter flow (L/min)	0.5 to 5 LPM
Sound level	<40dBA typical
Power Consumption (W)	<330 typical
Weight	18Kgs +/- 1Kg
Dimensions (mm)	410 x 265 x 566
Fuse	F5AL250V
Electrical Safety classification	Class II, BF type applied part; non AP/APG equipment
Ingress of water or particulate matter into equipment	IP21

Based on oxygen concentrations measure under STPD (standard temperature and pressure, 101.3kPa at an operating temperature of 20°C, dry).

Included with device

Filter	x4
Flat sponge air inlet filter	x1
Humidifier bottle	x1

Oxygen Concentration Flow Rate and Part Numbers

Oxygen Concentrator's flowrate	Oxygen concentration
0.5 L/min	90% - 96%
1 L/min	90% - 96%
2 L/min	90% - 96%
3 L/min	90% - 96%
4 L/min	90% - 96%
5 L/min	90% - 96%

Contact your local Philips Respironics sales representative regarding availability

Part number	Name of device
1156380	Oxygenate 5, KSW-5 230V 50Hz EU
1156381	Oxygenate 5, KSW-5 230V 50Hz UK
1156990	Oxygenate 5, KSW-5 230V 50Hz CH
1157200	Oxygenate 5, KSW-5 230V 60Hz SAUDI
1157201	Oxygenate 5, KSW-5 230V 50Hz UAE

Part number	Name of part
1156430	Oxygenate 5, Humidifier Bottle
1156429	Oxygenate 5, Filter Sponge (filter)
1157072	Oxygenate 5, Filter Sponge, Flat Foam (flat sponge air inlet filter)

* Air inlet filter, air filter, intake muffler filter, absorption tower intake filter, sieve filter, absorption tower effluent filter, oxygen tank filter, oxygen output filter and humidification bottle filter

Oxygen Concentrator – Oxygenate 5 is manufactured by Jiangsu Jonsung Bio-Medical Science and Technology Co. Ltd., distributed by Philips North America, LLC
© 2025 Koninklijke Philips N.V. All rights reserved.
Specifications are subject to change without notice.
www.philips.com
EM September 2025 MCI 4112764 | CAH 30/09/2025 |
00001255-08-v01 | Oxygenate 5, Datasheet, EMEA
1010 Murry Ridge Lane, Murrysville, PA 15668



Available in select markets.
Not available for sale in the United States.
Not all features available in all markets.